

## autoHSP - script environment

**WARNING:** You should **ONLY** use this script as a **script** to run. This script is not fit for running in the notebook environment. There will be abnormalities.

**NOTE:** When creating a new template/run script with this template, you should uncheck the `Pause on Warnings` options and change the timeout constraint to 5 hours.

```
If[True,
  Pause[1];
  If[
    $NotebookPage[TimeConstraint] < 4 Hour || Not[$NotebookPage[IgnoreWarnings]],
    PauseScript[];
    Return["You must increase the `TimeConstraint`
      to 5 hours and turn off `Pause on Warnings`."]
  ]
]
```

---

## Utility functions for logging

```
initVarGlobalLogMgr[reinit_ : False] := Block[{},
  If [! ValueQ[globalLogMgr] || reinit,
    globalLogMgr = {
      {"Output", "perm"}, {"Welcome to the HSP script developed by sijie."}
      (* tag, content *)
    };
  ];
];

initVarGlobalLogMgr[False];

getLogWithStyle[include_ : Nothing,
  exclude_ : Nothing, slice_ : All, returnPositionList_ : False,
  includeHow_ : SubsetQ, excludeHow_ : DisjointQ] := Block[
  {shouldHave, shouldNotHave, labels, logs, positions},
  shouldHave = If[ListQ[include], include, {include}];
  shouldNotHave = If[ListQ[exclude], exclude, {exclude}];
  positions = Range[Length[globalLogMgr]];
  labels = Quiet@Check[If[ListQ[#], #, {}] & /@ globalLogMgr[[All, 1]], {}];
  positions =
    Quiet@Check[Select[positions, includeHow[labels[[#]], shouldHave] &], {}];
  positions =
    Quiet@Check[Select[positions, excludeHow[labels[[#]], shouldNotHave] &], {}];
```

```

    positions = Quiet@Check[positions[[slice]], {}];
    logs = Quiet@Check[globalLogMgr[[positions]][All, 2], {}];
    logs = Quiet@Check[If[ListQ[#], StringJoin[ToString /@ #], If[StringQ[#],
        #, ToString[#]]], " --- Invalid input log --- "] & /@ logs;
    If[returnPositionList != True,
        Return[logs];,
        Return[positions];
    ];
];

printLogWithStyle[include_ : Nothing, exclude_ : Nothing, slice_ : All] := Block[{},
    Return[StringRiffle[getLogWithStyle[include, exclude, slice, False], "\n"]];
];

delLogWithStyle[include_ : Nothing, exclude_ : Nothing, slice_ : All] :=
    Block[{positions},
        positions = getLogWithStyle[include, exclude, slice, True];
        globalLogMgr = Quiet@Delete[globalLogMgr, List /@ positions];
    ];

InNotebookEnvQ[check_ : Cells] := Block[{},
    Return[Quiet@Check[ListQ[check[]], False]];
];

InScriptEnvQ[check_ : Cells] := Block[{},
    Return[NotebookPage[Type] === Object[Notebook, Script] &&
        Not[InNotebookEnvQ[check]]];
];

(* use the `smartCellPrint` function
to label your prints for future removals *)
smartCellPrint[objs_, labels_ : "tmp", styles_ : "Output", updateLog_ : True] :=
    Block[{content, cellTypes},
        content = StringJoin[ToString /@ If[ListQ[objs], objs, {objs}]];
        cellTypes = Join[ToString[#] & /@ If[ListQ[styles], styles, {styles}],
            ToString[#] & /@ If[ListQ[labels], labels, {labels}]];
        If[updateLog === True, AppendTo[globalLogMgr, {cellTypes, content}]];
        If[InNotebookEnvQ[], (* in a notebook environment *)
            CellPrint@Cell[BoxData[content], Sequence @@ cellTypes];,
            Return[content];
        ];
    ];

getCellWithStyle[include_ : Nothing, exclude_ : Nothing, slice_ : All] :=

```

```

Block[{shouldHave, shouldNotHave, cells},
  shouldHave = If[ListQ[include], include, {include}];
  shouldNotHave = If[ListQ[exclude], exclude, {exclude}];
  cells = Cells[FrontEnd`EvaluationNotebook[]];
  cells = Quiet@
    Check[Select[cells, SubsetQ[CurrentValue[#, CellStyle], shouldHave] &], {}];
  cells = Quiet@Check[Select[cells,
    DisjointQ[CurrentValue[#, CellStyle], shouldNotHave] &], {}];
  Return[Quiet@Check[cells[[slice]], {}]];
];

delCellWithStyle[include_ : Nothing, exclude_ : Nothing, slice_ : All] := Block[{},
  If[InNotebookEnvQ[], (* in a notebook environment *)
    NotebookDelete[getCellWithStyle[include, exclude, slice]];
  ];
  delLogWithStyle[include, exclude, slice];
];

smartCellPrint[
  "Welcome to the HSP script developed by sijie.", "perm", "Output", False];

```

Welcome to the HSP script developed by sijie.

## Functions to initialize variables

### ``globalSolvents``

``globalSolvents`` stores all the solvents used for this experiment with a mapping

```

initVarGlobalSolvents[reinit_ : False] := Block[{},
  If[reinit === True,
    smartCellPrint[
      {"Warning: You are force-reloading `initVarGlobalSolvents`. This
        could cause unexpected behaviors. Make
        sure you know what you are doing."},
    Nothing, {"Message", "MSG", "reinit-warning"}];
  ];
  If[! ValueQ[globalSolvents] || reinit,
    globalSolvents = <|
      "S1" → Model[Sample, "Methyl Ethyl Ketone"],
      "S2" → Model[Sample, "Hexanes"],
      "S3" → Model[Sample, "Toluene, Reagent Grade"],
      "S4" → Model[Sample, "Propylene Carbonate"],
      "S5" → Model[Sample, "Acetonitrile, HPLC Grade"],
    >|
  ];

```

```

    "S6" → Model[Sample, "Methanol"],
    "S7" → Model[Sample, "Diacetone Alcohol, 98%"],
    "S8" → Model[Sample, "Cyclohexanol, 99%"],
    "S9" → Model[Sample, "Dimethylformamide, Reagent Grade"],
    "S10" → Model[Sample, "n-Butyl Acetate, ACS Grade"],
    "S11" → Model[Sample, "Ethanol, Reagent Grade"],
    "S12" → Model[Sample, "1-Propanol"],
    "S13" → Model[Sample, "Cyclohexane"],
    "S14" → Model[Sample,
      "Propylene Glycol 1-Monomethyl Ether 2-Acetate (PGMEA)"],
    "S15" → Model[Sample, "Tetrahydrofuran, Anhydrous"],
    "S16" → Model[Sample, "Dichloromethane, Reagent Grade"],
    "S17" → Model[Sample, "N-Methyl-2-pyrrolidone, Reagent Grade"],
    "S18" → Model[Sample, "gamma Butyrolactone"],
    "S19" → Model[Sample, "1,4-Dioxane, anhydrous"],
    "S20" → Model[Sample, "Ethylene Glycol Monomethyl Ether Acetate"],
    (* NOT in inventory *)
    "S21" → Model[Sample, "Diethylene Glycol Monoethyl Ether Acetate"],
    "S22" → Model[Sample, "Diethyl ether"],
    "S23" → Model[Sample, "Glycerol"],
    "S24" → Model[Sample, "Formamide"],
    "S25" → Model[Sample, "Milli-Q water"]
  |>;
];
If[! ValueQ[globalSolventsUnits] || reinit,
  globalSolventsUnits = <|
    (* default unit for solvent is `Milliliter` *)
    "" → 1 Milliliter
  |>;
];
If[! ValueQ[getSolventUnit[]] || reinit,
  getSolventUnit[code_ : ""] :=
    Lookup[globalSolventsUnits, code, globalSolventsUnits[""]];
];
];

```

## `globalResins`

`globalResins` stores all the resins used for this experiment with a mapping

```

In[*]:= initVarGlobalResins[reinit_ : False] := Block[{},
  If[reinit === True,
    smartCellPrint[
      {"Warning: You are force-reloading `initVarGlobalResins`. This

```

```

        could cause unexpected behaviors. Make
        sure you know what you are doing."},
    Nothing, {"Message", "MSG", "reinit-warning"}]];
];
If[! ValueQ[globalResins] || reinit,
  globalResins = <|
    "R1" → Model[Sample, "Desmodur N 3300"], (* NOT in inventory *)
    "R2" → Model[Sample, "Desmodur N 3500"], (* NOT in inventory *)
    "R3" → Model[Sample, "Desmodur N 31100"], (* NOT in inventory *)
    "R4" → Model[Sample, "Desmodur N 3400"], (* NOT in inventory *)
    "R5" → Model[Sample, "Desmodur N 100"], (* NOT in inventory *)
    "R6" → Model[Sample, "Desmodur N 3900"], (* NOT in inventory *)
    "R7" → Model[Sample, "Desmodur Z 2589"], (* solid, NOT in inventory *)
    "R8" → Model[Sample, "Desmophen NH 1220"],
    "R9" → Model[Sample, "Desmophen NH 1420"], (* NOT in inventory *)
    "R10" → Model[Sample, "Desmophen NH 1520"],
    "R11" → Model[Sample, "Desmophen 1100"],
    "R12" → Model[Sample, "Desmophen 1200"],
    "R-2" → Model[Sample, "Poly(Ethylene Glycol) - 8000 MW"],
    (* solid flakes, NOT in inventory,
    optional resin to use if none of the above resins are available *)
    "R-1" → Model[Sample, "Polysorbate 20"] (* viscous fluid (250-450 mPa.s),
    default resin to use if none of the above resins are available *)
  |>;
];
If[! ValueQ[globalResinsUnits] || reinit,
  globalResinsUnits = <|
    (* By default, resins also use `Milliliter` as the unit *)
    "" → 1 Milliliter,
    "R7" → 1 Gram,
    (* Note that you can also change to `1.2Gram` for scaling *)
    "R-2" → 1 Gram
  |>;
];
If[! ValueQ[getResinUnit[]] || reinit,
  getResinUnit[code_ : ""] :=
    Lookup[globalResinsUnits, code, globalResinsUnits[""]];
];
];

```

## `globalVessel`

`globalVessel` stores information about the container that will be used. The following criteria should be met:

- A total of ~10 mL of liquid will be added to the vessel.
- The vessel should hold at least 15 mL of liquid. Preferably larger than 20 mL.
- The vessel should be compatible with the following shaker/mixer of choice.
- Clear glass walls (no texts, no marks, etc) and pure color caps (no horizontal textures).

```
initVarGlobalVessel[reinit_ : False, checkShaker_ : True] := Block[{},
  If[reinit === True,
    smartCellPrint[
      {"Warning: You are force-reloading `initVarGlobalVessel`. This
        could cause unexpected behaviors. Make
        sure you know what you are doing."},
      Nothing, {"Message", "MSG", "reinit-warning"}];
  ];

  If[! ValueQ[globalVessel] || reinit,
    (* if any of the variables not already defined or force `reinit`*)
    globalVessel = Model[Container, Vessel,
      "20mL Glass Scintillation Vial, With Cone-Shaped Cap Liner"];
    If[reinit && checkShaker && ValueQ[initVarGlobalShaker[]] &&
      ValueQ[globalShaker], initVarGlobalShaker[reinit, globalVessel];];
  ];
];
```

`globalShaker`, `globalShakeTime`, etc.

You need to make sure that the shaker chosen here is compatible with the previous `globalVessel`. A quick way to check this is by using `AnyTrue[MixDevices[globalVessel, 10 Milliliter], # == globalShaker &]`.

If you need to suppress the vessel-shaker check, pass `False` to `checkVessel`. Otherwise by default, `checkVessel` will try to find `globalVessel` and use it.

```
initVarGlobalShaker[reinit_ : False, checkVessel_ : True] := Block[{vessel},
  If[reinit === True,
    smartCellPrint[
      {"Warning: You are force-reloading `initVarGlobalShaker`. This
        could cause unexpected behaviors. Make
        sure you know what you are doing."},
      Nothing, {"Message", "MSG", "reinit-warning"}];
  ];

  If[! ValueQ[globalShaker] || reinit,
    (* if any of the variables not already defined or force `reinit`*)
    (* globalShaker=Model[Instrument,Vortex,"20 mm Bottle Vortex Genie"]; *)
    (* globalShaker=
      Model[Instrument,Shaker,"Burrell Scientific Wrist Action Shaker"]; *)
```

```

globalShaker = Model[Instrument, Roller, "Enviro-Genie"];

(* Shaker-vessel compatibility check *)
vessel = False;
(* we will not check the shaker-
vessel compatibility unless checkVessel is properly configured *)
If[checkVessel === True,
  If[AllTrue[{globalVessel}, ValueQ],
    vessel = globalVessel;,
    smartCellPrint[
      "`checkVessel` is set to True but `globalVessel` has not been set.",
      "initVarGlobalShaker"];
  ];
];
If[Quiet@Check[checkVessel[Type] === Model[Container, Vessel], False],
  vessel = checkVessel;,
  If[checkVessel != False && checkVessel != True,
    smartCellPrint["Wrong argument passed to `checkVessel`. Should be
      one of True|False|Model[Container,Vessel,...]",
      "initVarGlobalShaker"];
  ];
];

If[Not[vessel === False],
  smartCellPrint[
    "Checking for shaker-vessel compatibility... This might take
      some time...", "initVarGlobalShaker-incompatible"];
  If[AnyTrue[
    MixDevices[vessel, vessel[MaxVolume] / 2], # === globalShaker[Object] &],
    smartCellPrint[
      {"Success: shaker `" , globalShaker, "` is compatible with vessel `" ,
        vessel, "`."}, "initVarGlobalShaker"];,
    smartCellPrint[{"Warning: the shaker `" ,
      globalShaker, "` is not compatible with your vessel `" , vessel,
        "` . Please match these choices properly."}, "initVarGlobalShaker"];
  ];,
  smartCellPrint["Shaker-vessel compatibility test NOT performed.",
    "initVarGlobalShaker"];
];

delCellWithStyle["initVarGlobalShaker-incompatible"];
]; (* shaker initiation complete *)

If[! ValueQ[globalShakeTime] ||

```

```

    ! ValueQ[globalShakeTimeSolvent] || ! ValueQ[globalShakeRate] || reinit,
    globalShakeTime = 12 Hour; (* How long to mix the sample *)
    globalShakeTimeSolvent = 2 Hour;
    globalShakeRate = 15 RPM; (* This might not be applicable to the shaker *)
];
If[! ValueQ[globalSitTime] || reinit,
    globalSitTime = 24 Hour;
    (* how long to sit the sample for before taking the image *)
];
];

```

## `globalCamera`

```

In[*]:= initVarGlobalCamera[reinit_ : False] := Block[{},
    If[reinit === True,
        smartCellPrint[
            {"Warning: You are force-reloading `initVarGlobalCamera`. This
             could cause unexpected behaviors. Make
             sure you know what you are doing."},
            Nothing, {"Message", "MSG", "reinit-warning"}];
    ];

    If[! ValueQ[globalCamera] || ! ValueQ[globalCameraShowVessel] || reinit,
        (* if any of the variables not already defined or force `reinit` *)
        globalCamera =
            Model[Instrument, SampleImager, "Emerald DSLR Camera Imaging Station"];
        globalCameraShowVessel = True (* when imaging, image container or not *)
    ];
];

```

## `globalTaskMgr` and `globalLogMgr`

`globalTaskMgr` will serve as a collection of all tasks (“prep”, “image”, “EOE”) claimed by different threads running. For example:

```

globalTaskMgr = <|
    "thread1" -> <|
        "task" -> "prep", (* prepare samples *)
        "taskId" -> "xxxxxx",
        "samples" -> {
            <|"resin" -> "R1", "solvent" -> {"S1"}, "ramount" -> 5, "samount" -> {5}, "label" -> "R1S1"|>,
            <|"resin" -> "R1", "solvent" -> {"S7"}, "ramount" -> 5, "samount" -> {5}, "label" -> "R1S7"|>
        }
    |>,
    "thread2" -> <|
        "task" -> "image", (* image samples *)

```



```
"taskId" -> "xxxxxx",  
"samples" -> {  
  "R1S1", "R1S7"  
}  
|>  
|>
```

```

In[*]:= initVarGlobalTaskMgr[reinit_ : False, forsure_ : False] := Block[{},
  If[reinit == True,
    smartCellPrint[
      {"Warning: You are force-reloading `initVarGlobalTaskMgr`. This
        could cause unexpected behaviors. Make
        sure you know what you are doing."},
      Nothing, {"Message", "MSG", "reinit-warning"}];
  ];

  If[! ValueQ[globalTaskMgr] || And[reinit, forsure],
    (* if any of the variables not already defined or force `reinit` *)
    (* You really should NOT reinit `globalTaskMgr` after the script starts *)
    globalTaskMgr = <|
      (* `thread` → `task` *)
      |>;
  ];

  If[! ValueQ[getTaskForThread[]] || ! ValueQ[delTaskForThread[]] || reinit,
    getTaskForThread[thread_] := (
      If[! ValueQ[globalTaskMgr], globalTaskMgr = <| |>;];
      Return[Lookup[globalTaskMgr, thread, <| |>]];
    );
    delTaskForThread[thread_] := (
      If[! ValueQ[globalTaskMgr], globalTaskMgr = <| |>;];
      globalTaskMgr = KeyDrop[globalTaskMgr, thread];
      Return[globalTaskMgr];
    );
  ];

  If[! ValueQ[parseAssocList[]] || ! ValueQ[getPosOfLength[]] || reinit,
    (* if any of the variables not already defined or force `reinit` *)
    parseAssocList[key_, assocList_, fillMissWith_ : Nothing] := Cases[assocList,
      assoc_ /; AssociationQ[assoc] => Lookup[assoc, key, fillMissWith]];
    getPosOfLength[listOfLists_, N_ : 1,
      filterFunc_ : GreaterEqual, lengthFunc_ : Length] := Select[
      Range[Length[listOfLists]], filterFunc[lengthFunc[listOfLists[[#]]], N] &]
  ];
];

```

`globalNextUrl`, `globalUploadUrl`, and `parseJsonResponse`

Variables for server communication

`parseJsonResponse`: Function to parse JSON response from HTTP request

```

In[*]:= initVarGlobalUrl[reinit_ : False] := Block[{},
  If[reinit == True,
    smartCellPrint[
      {"Warning: You are force-reloading `initVarGlobalUrl`. This could cause
        unexpected behaviors. Make sure you know what you are
        doing."}, Nothing, {"Message", "MSG", "reinit-warning"}];
  ];

  If[! ValueQ[globalDomainName] || ! ValueQ[globalApiKey] || reinit,
    (* if any of the variables not already defined or force `reinit` *)
    globalDomainName = "https://XXXXX.XXXXX.cmu.edu/api/";
    globalApiKey = "XXXXXXXXXXXXXXXXXXXX"; (* API KEY *)
  ];

  If[! ValueQ[globalNextUrl] ||
    ! ValueQ[globalUploadUrl] || ! ValueQ[globalNotifyUrl] || reinit,
    globalNextUrl = URLBuild[{globalDomainName, globalApiKey, "next"}];
    globalUploadUrl = URLBuild[{globalDomainName, globalApiKey, "upload"}];
    globalNotifyUrl = URLBuild[{globalDomainName, globalApiKey, "notify"}];
  ];

  If[! ValueQ[globalMaxRuns] || ! ValueQ[globalWaitInterval] || reinit,
    globalMaxRuns = 10;
    (* try communicate with the server for a max of 10 times *)
    globalWaitInterval = 60;
    (* each failure will cause the program to pause before next try *)
  ];

  If[! ValueQ[parseJsonResponse[]] || reinit,
    parseJsonResponse[response_] := (
      ResourceFunction["ToAssociations"][
        ImportString[response["Body"], "JSON"]
      ];
  ];
];

```

## Other global variables

```

In[*]:= initVarGlobalOthers[reinit_ : False] := Block[{},
  If[reinit === True,
    smartCellPrint[
      {"Warning: You are force-reloading `initVarGlobalOthers`. This
        could cause unexpected behaviors. Make
        sure you know what you are doing."},
      Nothing, {"Message", "MSG", "reinit-warning"}];
  ];

  If[! ValueQ[globalExpSite] || reinit,
    globalExpSite = Object[Container, Site, "ECL-CMU"];
  ];

  If[! ValueQ[globalSampleNamePrefix] || reinit,
    globalSampleNamePrefix = "Washburn_HSP_";
  ];

  If[
    ! ValueQ[addSampleNamePrefix[]] || ! ValueQ[delSampleNamePrefix[]] || reinit,
    addSampleNamePrefix[string_, prefix_ : True] := Block[{tmpPrefix},
      tmpPrefix = If[StringQ[prefix], prefix, If[ValueQ[globalSampleNamePrefix],
        ToString[globalSampleNamePrefix], Return[string];
        ""]];
      If[StringLength[tmpPrefix] == 0, Return[string]];
      Return[If[! StringStartsQ[string, tmpPrefix | "id:"],
        StringJoin[tmpPrefix, string], string]];
    ];
    delSampleNamePrefix[string_, prefix_ : True] := Block[{tmpPrefix},
      tmpPrefix = If[StringQ[prefix], prefix, If[ValueQ[globalSampleNamePrefix],
        ToString[globalSampleNamePrefix], Return[string];
        ""]];
      If[StringLength[tmpPrefix] == 0, Return[string]];
      Return[If[StringStartsQ[string, tmpPrefix],
        StringReplace[string, StartOfString ~~ tmpPrefix -> ""], string]];
    ];
  ];
];

```

a wrapper to init all variables

```
In[*]:= initAllVars[reinit_ : False] := (
  initVarGlobalSolvents[reinit];
  initVarGlobalResins[reinit];
  initVarGlobalVessel[reinit, False];
  (* will NOT check if the vessel is compatible with the shaker *)
  initVarGlobalShaker[reinit, globalVessel];
  (* will also check if the vessel is compatible with the shaker *)
  initVarGlobalCamera[reinit];
  initVarGlobalTaskMgr[reinit, False];
  (* you should really NOT reinit this variable *)
  initVarGlobalUrl[reinit];
  initVarGlobalOthers[reinit];
);

(* initialize the variables *)
initAllVars[False];
(* delCellWithStyle[{"Message","MSG","reinit-warning"}]; *)
```

```
getLogWithStyle["initVarGlobalShaker"]
```

## Functions for server communications

### Function to pull the next experiments

Get the next task for a thread with defined abilities for the thread.

```
getNextTaskFromServer[thread_, ability_ : {"prep", "image"},
  extraOutputTag_ : Nothing, maxRuns_ : True, waitInterval_ : True] := Block[
  {tmpUrl, tmpUrlResponse, tmpCounter,
    tmpMaxRuns, tmpWaitInterval, tmpPrintCellStyle},
  initVarGlobalUrl[False];
  (* in case URL parameters are not initialized *)
  tmpMaxRuns = If[IntegerQ[maxRuns] && maxRuns ≥ 1, maxRuns, globalMaxRuns, 10];
  tmpWaitInterval = If[RealValuedNumberQ[waitInterval] && waitInterval > 0,
    waitInterval, globalWaitInterval, 60];
  tmpPrintCellStyle = {"getNextTaskFromServer", If[StringQ[extraOutputTag] &&
    StringLength[extraOutputTag] ≥ 1, extraOutputTag, Nothing, Nothing]};
  smartCellPrint[{"Trying to get the next task for thread `",
    thread, "` with abilities `", ability, "`"}, tmpPrintCellStyle];
  tmpUrl = URLBuild[globalNextUrl, <|"thread" → thread, "ability" → ability|>];
```

```

Do[
  tmpUrlResponse = URLRead[tmpUrl];
  If[tmpUrlResponse["StatusCode"] == 200,
    (
      (* If GET request successful, update task manager *)
      smartCellPrint[{"Try ", tmpCounter, "/", tmpMaxRuns,
        ": Trying to get the next task for thread `", thread,
        "` with abilities `", ability, "`"}, tmpPrintCellStyle];
      globalTaskMgr[thread] = jsonResponse[tmpUrlResponse];
      Break[];
    ),
    (* if not successful, pause 30 minutes until next retry *)
    smartCellPrint[{"Try ", tmpCounter, "/",
      maxRuns, " failed with code ", tmpUrlResponse["StatusCode"],
      ". Waiting ", tmpWaitInterval, " seconds before the next try."},
      {"getNextTaskFromServer", globalSmartCellPrintTag}];
    If[tmpCounter < tmpMaxRuns, Pause[tmpWaitInterval]];
  ],
  {tmpCounter, Max[1, tmpMaxRuns]}
]; (* Do ends *)

If[tmpUrlResponse["StatusCode"] != 200,
  smartCellPrint[{tmpMaxRuns, " URL GET request(s) failed with code ",
    tmpUrlResponse["StatusCode"], ": ", tmpUrl}, tmpPrintCellStyle];
  Return[False];
  delCellWithStyle[tmpPrintCellStyle];
  Return[True];
]
];

```

## Function to upload an image

For example, `uploadSampleImage["February 2024 User Training HPLC Sample", "test", 3, 1, False]`

```

uploadSampleImageToServer[sampleName_,
  taskId_, extraOutputTag_ : Nothing, validateTime_ : True,
  prefix_ : True, maxRuns_ : True, waitInterval_ : True] := Block[
  {prefixedSampleName, rawSampleName,
    tmpObjectSample, tmpObjectVessel, imageFromVessel, tmpSampleImage,
    tmpSampleImageTime, tmpSampleMixedN, tmpSampleMixedNTarget,
    tmpSampleMixedTime, tmpUrl, tmpPostData, tmpUrlResponse,
    tmpUrlResponseAsso, tmpMaxRuns, tmpWaitInterval, tmpPrintCellStyle},
  initVarGlobalUrl[False];
  (* in case URL parameters are not initialized*)

```

```

prefixedSampleName = addSampleNamePrefix[sampleName, prefix];
rawSampleName = delSampleNamePrefix[sampleName, prefix];
tmpMaxRuns = If[IntegerQ[maxRuns] && maxRuns ≥ 1, maxRuns, globalMaxRuns, 10];
tmpWaitInterval = If[RealValuedNumberQ[waitInterval] && waitInterval > 0,
  waitInterval, globalWaitInterval, 60];
tmpPrintCellStyle = {"uploadSampleImageToServer", If[StringQ[extraOutputTag] &&
  StringLength[extraOutputTag] ≥ 1, extraOutputTag, Nothing, Nothing]};

tmpObjectSample = Object[Sample, prefixedSampleName];
tmpObjectVessel = Object[Container, Vessel, prefixedSampleName];
(* issue in binding appearance log? *)
If[Length[tmpObjectVessel[AppearanceLog]] == 0 &&
  Length[tmpObjectSample[AppearanceLog]] == 0,
  smartCellPrint[{"Sample `", prefixedSampleName, "` from task `",
    taskId, "` does NOT yet have image data."}, tmpPrintCellStyle];
Return[False]; (* appearance data not ready yet *)
];
If[Length[tmpObjectVessel[AppearanceLog]] > 0 &&
  Length[tmpObjectSample[AppearanceLog]] > 0,
  imageFromVessel = If[tmpObjectVessel[AppearanceLog][[-1]][[1]] ≥
    tmpObjectSample[AppearanceLog][[-1]][[1]], True, False, False];,
  imageFromVessel =
    If[Length[tmpObjectSample[AppearanceLog]] == 0, True, False, True];
];

tmpSampleImage = If[imageFromVessel, tmpObjectVessel[AppearanceLog][[-1]][[2]][
  Image], tmpObjectSample[AppearanceLog][[-1]][[2]][Image]];
If[Head[tmpSampleImage] == Image,
  tmpSampleImageTime = If[imageFromVessel, tmpObjectVessel[AppearanceLog][[-1]][
    1], tmpObjectSample[AppearanceLog][[-1]][[1]],
  (* this should be the timestamp for the appearance *)
  smartCellPrint[{"Sample `", prefixedSampleName, "` from task `",
    taskId, "` does NOT yet have image data."}, tmpPrintCellStyle];
Return[False]; (* appearance data not ready yet *)
];

(* Image was successfully recorded,
next check if the image was taken after mixing *)
tmpSampleMixedN = 0;
tmpSampleMixedTime = tmpObjectSample[SampleHistory][[-1]][Date];

Do[
  If[Head[tmpOperation] == Mixed ||
    Quiet@Check[tmpOperation[Type] === Object[Protocol, Incubate], False],

```

```

    tmpSampleMixedN += 1;
    tmpSampleMixedTime = tmpOperation[DateCompleted];
  ], {tmpOperation, tmpObjectSample[Protocols]}
];
(* check if the sample was actually mixed *)
tmpSampleMixedNTarget = 1 +
  Boole[Length[StringCases[rawSampleName, RegularExpression["S\\d+"]] > 1];
If[tmpSampleMixedN < tmpSampleMixedNTarget && validateTime === True,
  smartCellPrint[{"Sample `", prefixedSampleName,
    "` from task `", taskId, "` has not been mixed properly for ",
    tmpSampleMixedNTarget, " times."}, tmpPrintCellStyle];
  Return[False];
];

(* check if the image time was after the mixed time *)
If[tmpSampleImageTime - tmpSampleMixedTime < globalSitTime &&
  validateTime === True,
  smartCellPrint[{"Sample `", prefixedSampleName, "` from task `",
    taskId, "` has not yet been imaged after a sit time of ",
    globalSitTime, " after mixing."}, tmpPrintCellStyle];
  Return[False];
];

(* All things checked, this is a valid sample image. Now upload it. *)
initVarGlobalUrl[False]; (* just in case *)
tmpUrl = URLBuild[globalUploadUrl, <|"sample" → rawSampleName,
  "sid" → tmpObjectSample[ID], "tid" → taskId|>];
tmpPostData = <|
  Method → "POST",
  "Body" → {
    "file" → <|
      "Content" → ExportByteArray[tmpSampleImage, "JPG"],
      "MIMEType" → "image/jpg",
      "Name" → ToString[rawSampleName] <> ".jpg"
    |>
  }
|>;

smartCellPrint[{"Image successfully retrieved ",
  If[validateTime === True, "and ", "but NOT "],
  "checked if its timestamp is valid. Trying to upload sample image of `",
  prefixedSampleName, "` under task id `", taskId,
  "` to the server."}, tmpPrintCellStyle];
Do[

```



```

smartCellPrint[{"Try ", tmpCounter, "/", tmpMaxRuns,
  ": Trying to upload sample image of `", prefixedSampleName,
  "` under task id `", taskId, "` to the server."}, tmpPrintCellStyle];
tmpUrlResponse = URLRead[HTTPRequest[tmpUrl, tmpPostData]];
If[tmpUrlResponse["StatusCode"] == 200,
  tmpUrlResponseAsso = parseJsonResponse[tmpUrlResponse];
  If[AssociationQ[tmpUrlResponseAsso] && Length[tmpUrlResponseAsso] ≥ 1 &&
    StringEndsQ[Keys[tmpUrlResponseAsso][[1]], ".jpg"],
    Break[];
  ];
];
(* if not successful, pause 30 minutes until next retry *)
smartCellPrint[{"Try ", tmpCounter, "/", tmpMaxRuns, " failed with code ",
  tmpUrlResponse["StatusCode"], ". Waiting ", tmpWaitInterval,
  " seconds before the next try."}, tmpPrintCellStyle];
If[tmpCounter < tmpMaxRuns, Pause[tmpWaitInterval];],
{tmpCounter, tmpMaxRuns}
];
If[tmpUrlResponse["StatusCode"] ≠ 200,
  smartCellPrint[{tmpMaxRuns, " URL POST request(s) failed with code ",
    tmpUrlResponse["StatusCode"], ": ", tmpUrl}, tmpPrintCellStyle];
  Return[False];,
  delCellWithStyle[tmpPrintCellStyle];
  Return[True];
]
];

```

Function to notify the server about experimentation status

```

notifyServer[params_, extraOutputTag_ : Nothing,
  maxRuns_Integer : True, waitInterval_ : True] := Block[
  {tmpUrl, tmpUrlResponse, tmpMaxRuns, tmpWaitInterval, tmpPrintCellStyle},
  initVarGlobalUrl[False]; (* just in case *)
  tmpUrl = If[AssociationQ[params], URLBuild[globalNotifyUrl, params], params];
  tmpMaxRuns = If[IntegerQ[maxRuns] && maxRuns ≥ 1, maxRuns, globalMaxRuns, 10];
  tmpWaitInterval = If[RealValuedNumberQ[waitInterval] && waitInterval > 0,
    waitInterval, globalWaitInterval, 60];
  tmpPrintCellStyle = {"notifyServer", If[StringQ[extraOutputTag] &&
    StringLength[extraOutputTag] ≥ 1, extraOutputTag, Nothing, Nothing]};

  Do[
    smartCellPrint[{"Try ", tmpCounter, "/", tmpMaxRuns,
      ": Trying to reach the server at: ", tmpUrl}, tmpPrintCellStyle];
    tmpUrlResponse = URLRead[tmpUrl];
    If[tmpUrlResponse["StatusCode"] == 200,
      Break[];
      (* if not successful, pause 30 minutes until next retry *)
      smartCellPrint[{"Try ", tmpCounter, "/", tmpMaxRuns, " failed with code ",
        tmpUrlResponse["StatusCode"], ". Waiting ", tmpWaitInterval,
        " seconds before the next try."}, tmpPrintCellStyle];
      If[tmpCounter < tmpMaxRuns, Pause[tmpWaitInterval]]];
    ];,
    {tmpCounter, tmpMaxRuns}
  ];

  delCellWithStyle[tmpPrintCellStyle];
  If[tmpUrlResponse["StatusCode"] ≠ 200,
    smartCellPrint[{tmpMaxRuns, " runs to reach the server failed with code ",
      tmpUrlResponse["StatusCode"], ": ", tmpUrl}, tmpPrintCellStyle];
    Return[False];,
    Return[True];
  ];
];

```

---

## Functions to run experiments

### Function to prepare sample solutions

The input to this function should be the Association coming from globalTaskMgr["thread\*\*\*"]

```

prepSamples[tasks_, extraOutputTag_ : Nothing, prefix_ : True] := Block[
  {falseReturnVal, intendedTask, resins,

```

```

    resinsAmounts, solvents, solventsAmounts, labels, taskList,
    prepProtocol, solventsIndices, validExpFlag, tmpPrintCellStyle},
(* You should have initiated some of
the global variables before you call this function *)
falseReturnVal = {False, False, False};
(* return this if an error occurs, {'status','protocol','samples'} *)
tmpPrintCellStyle = {"prepSamples", If[StringQ[extraOutputTag] &&
    StringLength[extraOutputTag] ≥ 1, extraOutputTag, Nothing, Nothing]};
smartCellPrint[
    {"Initiating function to build a prep sample protocol for tasks:", tasks},
    tmpPrintCellStyle];
intendedTask = Lookup[tasks, "task", ""];
If[! StringMatchQ[intendedTask, "prep", IgnoreCase → True],
    smartCellPrint[
        {"Wrong `prepSamples` function called to excute your intended task `",
        intendedTask, "`."}, tmpPrintCellStyle];
    Return[falseReturnVal];
];

(* preprocess requests in tasks *)
resins = parseAssocList["resin", tasks["samples"], False];
resinsAmounts = parseAssocList["ramount", tasks["samples"], False];
solvents = parseAssocList["solvent", tasks["samples"], False];
solventsAmounts = parseAssocList["samount", tasks["samples"], False];
labels = parseAssocList["label", tasks["samples"], False];
taskList = {resins, resinsAmounts, solvents, solventsAmounts, labels};
If[AnyTrue[taskList, MemberQ[#, False] &],
    smartCellPrint[
        {"Error: Your input tasks have missing entries in its members."},
        tmpPrintCellStyle];
    Return[falseReturnVal];
];

labels = addSampleNamePrefix[#, prefix] & /@ labels;
(* These labels will be `UploadName` names if not ids *)
If[! DuplicateFreeQ[labels],
    smartCellPrint[{"Error: Duplicated labels detected."}, tmpPrintCellStyle];
    Return[falseReturnVal];
];

If[Not[AllTrue[Flatten[resins], KeyExistsQ[globalResins, #] &]] ||
    Not[AllTrue[Flatten[solvents], KeyExistsQ[globalSolvents, #] &]],
    smartCellPrint[
        {"Error: Undefined resin(s) or solvent(s) found. Please check

```

```

        your inputs or update you `globalResins`
        and/or `globalSolvents`."}, tmpPrintCellStyle];
Return[falseReturnVal];
];

(* Start manual sample preparation protocol *)
smartCellPrint[{"A valid set of tasks"}, tmpPrintCellStyle];
Quiet@Check[
  prepProtocol = ManualSamplePreparation[
    (* 1. Label containers *)
    LabelContainer[
      Label → labels,
      Container → Table[globalVessel, Length[labels]]
    ],

    (* 2. Sequentially add the solvents into the vessels; note that some
       only have 1 solvent, while the others might have 2, 3, ... *)
    Sequence@@Reap[
      Do[
        solventsIndices =
          getPosOfLength[solvents, nSolvents, GreaterEqual, Length];
        If[
          Length[solventsIndices] < 1,
          Sow[Nothing]; Break[]; (* Upper bound reached *)
        ];
        Sow[
          Transfer[
            Source → Lookup[globalSolvents,
              solvents[solventsIndices][[All, nSolvents]]],
            Destination → labels[solventsIndices],
            Amount → solventsAmounts[solventsIndices][[All, nSolvents]] *
              getSolventUnit[solvents[solventsIndices][[All, nSolvents]]]
          ]
        ],
        {nSolvents, 10}
      ] (* using 10 solvents would be crazy if it ever happens*)
    ] [[2, 1]],

    (* 2+. Do we need to mix
       the samples with more than 1 solvent in them? *)
    solventsIndices = getPosOfLength[solvents, 2, GreaterEqual, Length];
    If[Length[solventsIndices] > 0,
      Mix[

```

```

    Sample → labels[solventsIndices],
    Instrument → globalShaker,
    MixUntilDissolved → False, (* Even though
        we'd want the solvents to dissolve with each other *)
    Time → globalShakeTimeSolvent,
    MixRate → globalShakeRate,
    ImageSample → True,
    MeasureVolume → False,
    MeasureWeight → False
],
Nothing (* Do nothing if all samples require at most 1 solvent *)
],

(* 3. Transfer resin after solvent *)
Transfer[
    Source → Lookup[globalResins, resins],
    Destination → labels,
    Amount → resinsAmounts * getResinUnit[resins]
],

(* 4. Shake the resin-solvent mixture *)
Mix[
    Sample → labels,
    Instrument → globalShaker,
    MixUntilDissolved → False, (* Dissolving is not guaranteed *)
    Time → globalShakeTime,
    MixRate → globalShakeRate
],

ImageSample → True,
MeasureVolume → False,
MeasureWeight → False
];,
smartCellPrint[
    {"Some unexpected error occurred when building the prep sample
        protocol. Please check again."}, tmpPrintCellStyle];
Return[falseReturnVal];
];

delCellWithStyle[tmpPrintCellStyle];
smartCellPrint[{"Starting to validate the generated prep sample protocol: ",
    prepProtocol}, tmpPrintCellStyle];
validExpFlag = Check[
    prepProtocol = ExperimentSamplePreparation[

```

```

        prepProtocol,
        ImageSample → True,
        MeasureVolume → False,
        MeasureWeight → False,
        OptimizeUnitOperations → True,
        Site → globalExpSite
    ];
    True,
    False
];
If[validExpFlag != True,
  smartCellPrint[
    {"Failed to generate a valid prep sample protocol for tasks: ", tasks},
    tmpPrintCellStyle];,
  delCellWithStyle[tmpPrintCellStyle, Nothing, -1];
];
Return[{validExpFlag, prepProtocol, labels}];
];

```

Function to upload sample names to the constellation from a prep sample protocol

```

rememberSamples[protocol_, samples_,
  taskId_ : "", extraOutputTag_ : Nothing, prefix_ : True] := Block[
  {whichProtocol, returnVal, names, samplesObjects, i, tmpPrintCellStyle},
  tmpPrintCellStyle = {"rememberSamples", If[StringQ[extraOutputTag] &&
    StringLength[extraOutputTag] ≥ 1, extraOutputTag, Nothing, Nothing]};
  whichProtocol = If[StringQ[protocol],
    Object[Protocol, ManualSamplePreparation, protocol], protocol];
  If[Head[taskId] === String && StringLength[taskId] > 0,
    Quiet@UploadName[whichProtocol, addSampleNamePrefix[taskId]];
  ];

  names = If[ListQ[samples], samples, {samples}];
  names = addSampleNamePrefix[#, prefix] & /@ names;
  (* These labels should be `UploadName` names if not ids *)

  smartCellPrint[{"Got inputs protocol `" , protocol, "` and samples `" , names,
    "` . Now trying to upload the sample names..."}, tmpPrintCellStyle];
  samplesObjects = Quiet@LookupLabeledObject[whichProtocol, names];
  returnVal = # === Null & /@ samplesObjects;
  For[i = 1, i ≤ Length[names], i++,
    returnVal[[i]] = Check[ListQ[UploadName[samplesObjects[[i]], names[[i]]], False];
  ];
];

```

```

(* Print some error messages *)
For[i = 1, i ≤ Length[names], i++,
  If[returnVal[[i]] == False,
    If[samplesObjects[[i]] == Null,
      smartCellPrint[{"Error: Could not find name `" , names[[i]],
        "` in protocol `" , protocol, "`."}, tmpPrintCellStyle];
      smartCellPrint[{"Error: Could not upload name `" , names[[i]],
        "` for the found sample object `" , samplesObjects[[i]],
        "` in protocol `" , protocol, "`."}, tmpPrintCellStyle];
    ];
    (* the sample name was uploaded successfully →
    Now upload its sample content and source protocol *)
    If[samplesObjects[[i]][Type] === Object[Sample],
      UploadName[samplesObjects[[i]][Container], names[[i]]];
    ];
    If[samplesObjects[[i]][Type] === Object[Container, Vessel],
      UploadName[samplesObjects[[i]][ContentsLog][[-1]][[3]], names[[i]]];
    ];
  ];
];

If[AllTrue[returnVal, TrueQ], delCellWithStyle[tmpPrintCellStyle];];
Return[returnVal];
];

```

## Function to image samples

```

imageSamples[tasks_, extraOutputTag_ : Nothing, prefix_ : True] := Block[
  {falseReturnVal, intendedTask,
    samples, imageProtocol, validExpFlag, tmpPrintCellStyle},
  falseReturnVal = {False, False, False};
  (* return this if an error occurs, {'status', 'protocol', 'samples'} *)
  tmpPrintCellStyle = {"imageSamples", If[StringQ[extraOutputTag] &&
    StringLength[extraOutputTag] ≥ 1, extraOutputTag, Nothing, Nothing]};
  smartCellPrint[
    {"Initiating function to build a image sample protocol for tasks:", tasks},
    "imageSamples"];

  intendedTask = Lookup[tasks, "task", ""];
  If[! StringMatchQ[intendedTask, "image", IgnoreCase → True],
    smartCellPrint[
      {"Wrong `imageSamples` function called to excute your intended task `",

```

```

        intendedTask, "`."), tmpPrintCellStyle];
    Return[falseReturnVal];
];

samples = Lookup[tasks, "samples", {}];
samples = If[! ListQ[samples], {samples}, samples];
samples = addSampleNamePrefix[#, prefix] & /@ samples;
(* These labels should be `UploadName` names if not ids *)
If[Not[DuplicateFreeQ[samples]],
    smartCellPrint[{"Found duplicates in your parsed sample names: `",
        samples, "` . Please check again."}, tmpPrintCellStyle];
    Return[falseReturnVal];
];

(* Next, search for these samples and image them *)
smartCellPrint[{"Starting to validate the generated image sample protocol: ",
    imageProtocol}, tmpPrintCellStyle];
validExpFlag = Check[
    ValidExperimentImageSampleQ[
        Object[Sample, #] & /@ samples,
        Instrument → globalCamera,
        ImageContainer → globalCameraShowVessel,
        Site → globalExpSite
    ] == True,
    smartCellPrint[
        {"Some unexpected error occurred when building the image sample
            protocol. Probabaly because one or more of the
            sample names are not correctly recognized. Please
            check again.", samples}, tmpPrintCellStyle];
    Return[falseReturnVal];
    False
];

If[validExpFlag != True,
    smartCellPrint[
        {"Failed to generate a valid image sample protocol for tasks: ", tasks},
        tmpPrintCellStyle];
    Return[falseReturnVal];
];

delCellWithStyle[tmpPrintCellStyle];
imageProtocol = ExperimentImageSample[
    Object[Container, Vessel, #] & /@ samples,
    Instrument → globalCamera,
    ImageContainer → globalCameraShowVessel,

```



```
Site → globalExpSite  
];  
  
Return[{validExpFlag, imageProtocol, samples}];  
];
```

---

## HSP - Main (putting things together)

Some preparatory work

```

thread = "thread1";
(* should not contain characters such as `-'`, `@` and `:` *)
ability = {"prep", "image"};
validateImageTime = True;
extraOutputTag = "HSPMain:" <> thread;
prefix = globalSampleNamePrefix; (* Washburn_HSP_*)
maxRuns = globalMaxRuns;
waitInterval = globalWaitInterval;

whichCycle = 0;
useLastTask = False;
(* if last cycle's task went wrong and you need to rerun it *)
thisCycleWentFine = True;
(* if everything is executing without error/exception *)
EOE = False; (* end of experiment *)

(* pauseExecution will handle errors and exceptions *)
pauseExecution[
  params_ : <|"status" → "exception",
    "message" → "Your attention is required ASAP."|>,
  extraOutputTag_ : Nothing, maxRuns_ : True, waitInterval_ : True] := (
  smartCellPrint[
    {"Pausing script execution due to error or exception.\nYou are
      advised to review the error information carefully and make
      adjustments on the server before resuming execution."},
    {"pauseExecution", extraOutputTag}];
  thisCycleWentFine = False;
  notifyServer[params, extraOutputTag, maxRuns, waitInterval];
  PauseScript[];
  );

printLineBreak[cycle_ : Nothing, extraOutputTag_ : Nothing] := (
  Return[smartCellPrint[{StringRepeat["-", 12], " ", cycle, " ",
    StringRepeat["-", 12]}, Flatten[{"perm", "line", extraOutputTag}]]];
  );

printLineBreak["START", extraOutputTag];
smartCellPrint[{"Start HSPMain (script) execution on thread `" , thread, "`."},
  {"perm", extraOutputTag}];

```

## Start cycle from here

```
In[*]:= Label["start:HSP"];

In[*]:= whichCycle += 1;
tasks = <|>; (* the set of tasks to be performed *)
taskType = ""; (* which task will be performed in this cycle *)
extraOutputTag = "HSPMain:" <> thread <> "@" <> ToString[whichCycle];
thisCycleWentFine = True; (* a fresh new start,
you'd better behave well this time *)
thisTaskGotMatched = False;
failMessage = <|"status" → "exception",
  "message" → "An exception happened. Your attention is required ASAP."|>;
printLineBreak[whichCycle, extraOutputTag];
```

### 1 - Check if there is already thread-claimed tasks in `globalTaskMgr`

this should rarely happen since before each `Goto`, the tasks are removed

```
If[thisCycleWentFine && And[KeyExistsQ[globalTaskMgr, thread], Not[useLastTask]],
  (* Why is this thread NOT removed from task manager? *)
  smartCellPrint[{"You are trying to run thread `", thread,
    "` but it is already registered under `globalTaskMgr` with tasks: ",
    getTaskForThread[thread],
    ".\nDid you forget to `delTaskForThread`? You should also check if the task
    is already successfully performed in the lab.\nWARNING: make sure to
    update the server configuration before you skip the pause."},
    {"threadExistsException", extraOutputTag}];
thisCycleWentFine = False;
failMessage = <|"status" → "exception", "message" → "threadExistsException"|>;
];
```

## 2 - Get next set of tasks from the server

```

If[thisCycleWentFine,
  If[And[useLastTask, KeyExistsQ[globalTaskMgr, thread]],
    taskStatus = True,
    taskStatus = getNextTaskFromServer[thread,
      ability, extraOutputTag, maxRuns, waitInterval] == True
  ];
  If[Not[taskStatus],
    thisCycleWentFine = False;
    failMessage = <|"status" → "exception",
      "message" → "getNextTaskFromServerException"|>;
    Return[smartCellPrint[{"Cycle ", whichCycle,
      ": failed to get the next set of tasks from the server."},
      {"getNextTaskFromServer", extraOutputTag}]]
  ,
  tasks = getTaskForThread[thread];
  taskType = Lookup[tasks, "task", ""];
  If[! StringMatchQ[taskType, "pause", IgnoreCase → True],
    Return[smartCellPrint[
      {"Cycle ", whichCycle, ": got the next set of tasks from the server: ",
        tasks}, {"getNextTaskFromServer", extraOutputTag}]]
  ]
]
]

```

### 3.1 - prepSamples (?)

```

If[thisCycleWentFine && StringMatchQ[taskType, "prep", IgnoreCase → True],
  thisTaskGotMatched = True;
  urlParams = <|"status" → "exception", "message" → "unknownException"|>;

  {status, protocol, samples} = prepSamples[tasks, extraOutputTag];
  If[status != True,
    thisCycleWentFine = False;
    failMessage = <|"status" → "exception", "message" → "prepSamplesException"|>;
  ]
]

```

Supposedly, if the previous cell returns a valid protocol, it should automatically run. And the next cell will only start executing when that protocol has a status of completed.

**NOTE:** If the above cell (prepSamples) runs successfully, but the following cell (rememberSamples) failed. You should go to a notebook page and fix `rememberSamples` from there. DO NOT try to fix it in the script environment.

```

If[thisCycleWentFine && status && StringMatchQ[taskType,
  "prep", IgnoreCase → True] && protocol[Status] === Completed,
  status = rememberSamples[protocol,
    samples, Lookup[tasks, "taskId", ""], extraOutputTag];

If[AllTrue[status, TrueQ],
  protocolProgress = protocol[CheckpointProgress];
  samplePrepEndTime = "";
  If[Quiet@ListQ[protocolProgress],
    Do[
      If[progress[[1]] === "Sample Preparation",
        samplePrepEndTime = DateString[progress[[3]], "ISODateTime"] <>
          DateString[progress[[3]], "ISOTimeZone"];
        Break[];
      ];
      , {progress, Reverse[protocolProgress]}
    ];
  ];

urlParams = <|"status" → "success", "task" → "prep",
  "taskId" → Lookup[tasks, "taskId", ""], "protocolId" → protocol[ID],
  "samples" → (delSampleNamePrefix[#, prefix] & /@ samples),
  "sampleIds" → (Object[Container, Vessel, #][ID] & /@ samples),
  "time" → samplePrepEndTime|>;
If[notifyServer[urlParams, extraOutputTag, maxRuns, waitInterval] != True,
  thisCycleWentFine = False;
  failMessage = <|"status" → "exception",
    "message" → "rememberSamplesNotifyServerException"|>;
];
,
thisCycleWentFine = False;
failMessage = <|"status" → "exception",
  "message" → "rememberSamplesException"|>;
];
];

```

### 3.2 - imageSamples (?)

```

If[thisCycleWentFine && StringMatchQ[taskType, "image", IgnoreCase → True],
  thisTaskGotMatched = True;
  urlParams = <|"status" → "exception", "message" → "unknownException"|>;

  {status, protocol, samples} = imageSamples[tasks, extraOutputTag];
  If[status != True,
    thisCycleWentFine = False;
    failMessage = <|"status" → "exception", "message" → "imageSamplesException"|>;
  ]
]

```

Supposedly, if the previous cell returns a valid protocol, it should automatically run. And the next cell will only start executing when that protocol has a status of completed.

**NOTE:** If the above cell (imageSamples) runs successfully, but the following cell (uploadSampleImageToServer) failed. You should go to a notebook page and fix uploadSampleImageToServer from there. DO NOT try to fix it in the script environment.

```

If[thisCycleWentFine && status && StringMatchQ[taskType,
  "image", IgnoreCase → True] && protocol[Status] === Completed,
  status =
    uploadSampleImageToServer[#, Lookup[tasks, "taskId", ""], extraOutputTag,
      validateImageTime, prefix, maxRuns, waitInterval] & /@ samples;
  If[AllTrue[status, TrueQ],
    urlParams = <|"status" → "success", "task" → "image",
      "taskId" → Lookup[tasks, "taskId", ""], "protocolId" → protocol[ID],
      "samples" → (delSampleNamePrefix[#, prefix] & /@ samples),
      "sampleIds" → (Object[Container, Vessel, #][ID] & /@ samples) |>;
    If[notifyServer[urlParams, extraOutputTag, maxRuns, waitInterval] != True,
      thisCycleWentFine = False;
      failMessage = <|"status" → "exception",
        "message" → "uploadSampleImageToServerNotifyServerException"|>;
    ];
    ,
    thisCycleWentFine = False;
    failMessage = <|"status" → "exception",
      "message" → "uploadSampleImageToServerException"|>;
    ];
  ];
]

```

### 3.3 - EOE (?)

```
If[thisCycleWentFine && StringMatchQ[taskType, "EOE", IgnoreCase → True],
  EOE = True;
  thisTaskGotMatched = True;
  Return[smartCellPrint[{"EOE string received at cycle ", whichCycle, "."},
    {"perm", "EOE", extraOutputTag}]]
]
```

### 3.4 - pause (?)

```
(* initialize some temporary variables *)
If[thisCycleWentFine && StringMatchQ[taskType, "pause", IgnoreCase → True],
  thisTaskGotMatched = True;
  whichCycle -= 1;
  donotuseme = 1;
  For[i = 1, i ≤ 3 * 60, i++,
    (* donotuseme=Mod[donotuseme*i+i,2024]; *)
    donotuseme =
      parseJsonResponse[URLRead[URLBuild[{globalDomainName, "pause"}]]];
  ];
  If[MemberQ[globalLogMgr[[-1, 1]], "line"],
    delCellWithStyle[{"perm", "line"}, Nothing, -1];
  ];
  If[! MemberQ[globalLogMgr[[-1, 1]], "pause"],
    printLineBreak["SCRIPT PAUSED", {"pause", extraOutputTag}];
  ];
];
```

### 3.5 - All other cases: thisTaskDidNotGotMatched?

what is going on on the server???

```
If[thisCycleWentFine && thisTaskGotMatched != True,
  thisCycleWentFine = False;
  failMessage = <|"status" → "exception", "message" → "unknownTaskException"|>;
];
```

## Exceptions? Errors? End of experiment?

```

If[thisCycleWentFine != True,
  pauseExecution[failMessage, extraOutputTag, maxRuns, waitInterval];
  smartCellPrint[{"Cycle ", whichCycle,
    " exited with an exception/error. Good luck fixing it!"},
    {"perm", extraOutputTag}];
Return[printLogWithStyle[extraOutputTag]];
If[! StringMatchQ[taskType, "pause", IgnoreCase → True],
  smartCellPrint[{"Cycle ", whichCycle, " exited with success. Congrats!"},
    {"perm", extraOutputTag}];
Return[printLogWithStyle[extraOutputTag]];
]
]

```

If you need to perform some task before going back to next loop, edit the following cell before it is too late!

```

(* your custom correction code starts here *)
useLastTask = False;
(* if you'd like to keep the last task, change to `True` *)

(* your custom correction code ends here *)

(* builtin error message cleaning code starts here*)
If[thisCycleWentFine != True,
  delCellWithStyle[{"threadExistsException", extraOutputTag}];
  delCellWithStyle[{"pauseExecution", extraOutputTag}];
  delCellWithStyle[{"getNextTaskFromServer", extraOutputTag}];
  delCellWithStyle[{"uploadSampleImageToServer", extraOutputTag}];
  delCellWithStyle[{"notifyServer", extraOutputTag}];
  delCellWithStyle[{"prepSamples", extraOutputTag}];
  delCellWithStyle[{"rememberSamples", extraOutputTag}];
  delCellWithStyle[{"imageSamples", extraOutputTag}];
];

(* delete the thread claimed task *)
If[Not[useLastTask],
  globalTaskMgr = delTaskForThread[thread];
];

```



Not end of experiment (EOE)? Go back to the start!

```
If[EOE != True,
  Goto["start:HSP"];,
  PauseScript[]; (* pause for final human confirmation *)
];
```

**EOE at last!** Let's print out the whole log.

```
extraOutputTag = "HSPMain:" <> thread;
printLineBreak["END", extraOutputTag];
smartCellPrint[{"End HSPMain (script) execution on thread `", thread,
  "`.\nPlease restart execution if needed be.\nThank you for using our
  automated HSP measurment script!"}, {"perm", extraOutputTag}];
```

```
log = printLogWithStyle[]
```